


One Click® Water Determination

Simple and Secure




One Click® User Interface

The large color touchscreen features the simple-to-understand and fast-to-operate One Click user interface. All users have their own homescreen with dedicated shortcut buttons to directly start specific routine tasks.




Solvent Manager

Users are protected against contact with chemicals using the Solvent Manager liquid handling device. It is fully controlled by the titrator and offers solvent exchange with only One Click without additional user interaction.



Intelligent Burettes

The titrant is automatically recognized, when the burette is installed on the titrator. Important data such as concentration or lifetime is saved on the RFID chip and directly transferred to the titrator, where it is stored and can be monitored. Using the wrong or an expired titrant is therefore impossible.



PC software LabX® Titration

Full control over the titrator, secure database archiving, networking, efficient result management with control charts and individual search filters or rapid and clear method development are only some aspects of the powerful LabX titration PC software.



V20S Compact Volumetric KF Titrator

The routine Karl Fischer titrator combines fast and precise water content determination with easy operation. Samples from any industry with water content between 100 ppm and 100% are analyzed in single measurements or sample series with statistics.

- Color touchscreen with One Click user interface
- User interface in English, German, French, Spanish, Italian, Chinese, Russian, Polish, Korean and Portuguese
- Intelligent burettes with RFID chip saving all relevant titrant data
- Connectivity to LabX Titration PC software via USB or direct Ethernet (TCP/IP)
- Common USB printers or the USB-P25 compact printer for secure and GLP-compliant paper recording

| Functionality Overview | | V20S |
|------------------------|---|------|
| Application | Solid, liquid and gaseous samples | • |
| | Resources and method database | • |
| | KF titrants | 100 |
| | KF standards | 50 |
| One Click® Titration | User specific Home Screen with Shortcuts | • |
| | Shortcuts per user | 12 |
| | Solvent Control | - |
| Hot Plug & Play | Solvent Manager | • |
| | Burette recognition with titrant, concentration and monitoring data | • |
| | External Dosing Unit for addition of auxiliary reagents | - |
| | USB printer DIN-A4/US letter / USB-P25 compact printer | • |
| | USB memory stick / USB barcode reader | • |
| Security | LevelSens | - |
| | LogStraight Fingerprint Reader | - |
| Automation | InMotion KF Oven Autosamplers | - |
| User Management | "Expert – Routine" user groups | • |
| | Preprogrammed METTLER TOLEDO methods | - |
| | KF standard methods | 4 |
| | Max. number of user methods | 5 |
| | External extraction/dissolution dedicated method type | - |
| | Number of method function | 12 |
| | Number of samples per series | 120 |
| | Autostart: automatic start after sample addition | • |
| | Drift stop for endpoint recognition (absolute and relative, tmax, delay) | • |
| | Max. number of calculations/results per method | 3 |
| | User defined calculations | - |
| | Statistics of up to 120 samples including result limits check, outlier test | • |
| | Extra statistical functions as max. srel and multiple determinations | - |
| | Recalculation of results | • |
| | Data in Online Screen: time, water count, potential, stirrer speed, volume, drift | • |
| | Data in measured value table: time, water count, potential, volume, drift | • |
| | Curves in Online Screen and printout: E-V, E-t, V-t, drift-t, water count-t | • |
| Task list | Number of tasks queued | - |
| Online Help | Context sensitive parameter explanation | • |
| Data Export | USB, Ethernet, PDF | • |
| Languages | User interface and printouts in English, German, French, Spanish, Italian, Chinese, Russian, Polish, Korean, Portuguese | • |
| PC software | LabX® Titration Express and Server | • |
| Homogenizer | TTL | - |
| | Ethernet (TCP/IP) | 1 |
| | COM (serial balance connection) | 1 |
| | USB host (printer, memory stick, barcode reader, hub) | 2 |
| | USB client (PC) | 1 |
| | Integrated Magnetic Stirrer | 1 |
| | Pump / Stirrer | 1 |

- exclusive V30S functionality

Technical Specifications

| | | |
|--|---|-------------------------------------|
| Polarized sensor input | Measurement range | -2000...+2000 mV |
| | Resolution / Limits of error | 0.1 mV / 2 mV |
| | Polarization current range | 0...24 µA |
| | Polarization current resolution / Limits of error | 0.1 µA / 1 µA |
| Burette drive | Drive resolution | 20 000 steps |
| | Burette resolution (for 5 mL burette) | 0.25 µL |
| | Limits of error | 0.3% of the respective volume |
| | Resolver | 32 pulses / mm stroke |
| | Resolver resolution | 0.0625% of burette volume |
| | Filling time and discharge time | 20 s at 100% filling rate |
| | | |
| Environment | Ambient temperature | +5...+40 °C |
| | Atmospheric humidity | max. 80% at 31 °C |
| Titrator dimensions | Width x depth x height / Weight | 210 x 333 x 320 mm / 4.2 kg |
| | Glass beaker | 150 mL |
| Control and display | Touchscreen | VGA 5.7" color TFT, 640 x 480 pixel |
| Application (depends on the environment and handling) | Drift measurement | Online, < 5 µg/min |
| | Measurement range | 100 ppm...100% |
| | Repeatability | 0.3% at >10 mg H ₂ O |

METTLER TOLEDO Group

Laboratory Division

Local contact: www.mt.com/contact

Subject to technical changes

© 06/2017 METTLER TOLEDO. All rights reserved.

30265532B

Marketing Titration / MarCom Analytical

www.mt.com/titration

For more information



One Click is a Registered Trademark of METTLER TOLEDO in Switzerland, the European Union, Russia and Singapore.
LabX is a Registered Trademark of METTLER TOLEDO in Switzerland, USA, China, Germany and a further 13 countries.