

500 μL, Model 1750 RN SYR, Large Removable NDL, 22 ga, 2 in, point style 2

Specification Sheet (Part/REF # 81230)



Split/splitless injection consists of a heated chamber with a glass liner into which the sample is injected through the septum. A microsyringe is used to inject the sample through a rubber septum into a flash vaporizer chamber at the head of the column.

Product Specifications

Volume	500 μL
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Large RN
Gauge	22 gauge
Point Style	2
Needle Length	2 inch (51 mm)
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	700 psig (48.3 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.128 inches (3.26 mm)



Visit <u>www.hamiltoncompany.com/support</u> to find a representative in your area or call us:

US: +1-775-525-3000 EU: +40-356-635-055



Barrel Outer Diameter	0.305 inches (7.75 mm)
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std Needle Dead Volume	6.81 μL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice https://www.hamiltoncompany.com/laboratory-products/syringes/81230

Spec. Version B



Visit www.hamiltoncompany.com/support to find a representative in your area or call us:

US: +1-775-525-3000 EU: +40-356-635-055